



CLISIS *Evolution*

R/F Remote-controlled
system with digital flat
panel detector

THE EVOLUTION OF A TOP QUALITY PRODUCT

REDESIGNED, REENGINEERED, RESTYLED TO OFFER YOU
THE BEST EXPERIENCE EVER

OUTSTANDING PERFORMANCES

CLISIS is the **evolution** of a remote-controlled table, representing the best combination of operational **functionality**, premium **design** and enhanced image **quality** in any radiographic and fluoroscopic procedure

- **Advanced UI**
- **New design**
- Outstanding operational **efficiency**
- Complete **reliability**
- State-of-the-art **digital flat panel detector**
- Powerful **image processing**
- Fully **DICOM** compliant





TOTAL ACCESSIBILITY

The rear accessibility offered by the CLISIS tabletop allows easy patient positioning, granting maximum **comfort** and **efficiency**.



ALL IN ONE:
FUNCTIONALITY AND DESIGN





- **90/90** remote tilting system
- Wide adjustable **SID**
- **Reduced** tabletop-floor **distance**
- Single-end suspended carbon fiber tabletop
- Enhanced maximum patient **weight**
- Wide displacement of tube/ detector assembly for a **full-length patient coverage**



WIDE RANGE OF APPLICATIONS

The new CLISIS brings to **higher innovation** and **optimized workflow** both in routine and specialized procedures: continuous and pulsed fluoroscopy, high resolution digital spot exposures, tomography, tomosynthesis, spine and lower limbs automatic image stitching.



OBLIQUE PROJECTION



EXAMINATIONS IN
TRENDELENBURG POSITION



UROGENITAL STUDIES



SPECIAL PROJECTIONS
ON STRETCHER



SPECIAL PROJECTIONS WITH
WIRELESS DETECTOR



GASTROINTESTINAL
EXAMINATIONS



WIDE RANGE OF ACCESSORIES

UTMOST VERSATILITY FOR SPECIALIZED APPLICATIONS

Easy examinations of patients
with **reduced mobility**



EXAMINATIONS WITH HEIGHT ADJUSTABLE FOOTREST



EMERGENCY APPLICATIONS: SAFETY AND AGILITY

A smooth stretchers' transfer ensures **immediate intervention** also in emergency applications.



INNOVATIVE USER INTERFACE

New integrated console:
complete remote and quick control of the overall system



**Redesigned ergonomics
and intuitive user interface:**
maximum efficiency for the operator





DIAGNOSIS: QUALITY AND ACCURACY

Significantly improved processes, dose reduction and image enhancement is the result of a **renovated workflow** from the image acquisition to its storage and transmission.

Each specific exam with fluoro, rad, tomography and tomosynthesis procedure automatically elaborates processing algorithms, reducing **waiting times** and displaying **optimized images** almost in real time

MINIMUM PATIENT DOSE

Ideally suited for **pediatric procedures** thanks to the minimum dose delivered and the possibility to work in **direct contact** with the detector



GET IN TOUCH

GENERAL MEDICAL MERATE SPA
Via Partigiani, 25
24068 SERIATE (BG) - Italy

Phone +39 035 45 25 311
Fax +39 035 29 77 87
info@gmmspa.com



GMM
GROUP

www.gmmspa.com